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ایمیل: [info@novingene.com](mailto:info@novingene.com)

وبسایت: [www.novingene.com](http://www.novingene.com)

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- Münz, C. ed., 2015. Epstein Barr virus volume 2: one herpes virus: many diseases (Vol. 391). Springer.
- Kimura, H. and Kwong, Y.L., 2019. EBV viral loads in diagnosis, monitoring, and response assessment. *Frontiers in oncology*, 9, p.62.
- Yin, H., Qu, J., Peng, Q. and Gan, R., 2019. Molecular mechanisms of EBV-driven cell cycle progression and oncogenesis. *Medical microbiology and immunology*, 208, pp.573-583.
- Mackay, Ian M. "Real-time PCR in the microbiology laboratory." *Clinical microbiology and infection* 10, no. 3 (2004): 190-212.

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# EBV RQ Kit Manual

Summer 2024, Version 1.1

For Real-Time PCR Detection and Quantitation of Epstein-Barr virus (EBV) DNA

For use with Rotor-Gene or StepOne

For Research Use Only

 24 (Cat# EBVRQ24)

 48 (Cat# EBVRQ48)

 96 (Cat# EBVRQ96)


 NG-WI-ASL-32-101

RUO

**NovinGene ParsVira**

No. 5, Iranshahr St, Tehran, Iran 1581633336.

# Table of Contents

1. Introduction .....	3
2. Intended Use .....	3
3. Background Information .....	3
4. Test Principle .....	4
5. Kit Contents .....	4
6. Packaging models .....	4
7. Storage and Stability .....	5
8. Product Use Limitations .....	5
9. Additionally Required Materials .....	5
10. General Precautions .....	6
11. Specimen, Storage and Transport .....	6
12. Interfering Substances .....	6
13. Internal Control (IC) .....	7
14. DNA Isolation .....	7
15. PCR Protocol .....	8
16. Devices and software .....	8
17. Programming Rotor-Gene .....	8
18. Programming StepOne .....	9
19. Programming Other Machines .....	9

20. Data Analysis: Rotor-Gene .....	9
21. Data Analysis: StepOne .....	12
22. Quantitation .....	14
23. Linear Range .....	14
24. Sensitivity.....	14
25. Disposal Method .....	15
26. Technical Support.....	15
27. Contact Information.....	15
28. References .....	15
29. Symbols.....	16

## 1. Introduction

EBV RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting and quantifying Epstein-Barr virus (EBV) DNA. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

## 2. Intended Use

EBV RQ kit is intended for detecting and quantifying Epstein-Barr virus (EBV) DNA. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene or StepOne machines.

## 3. Background Information

Human Herpesvirus 4 (HHV-4) known as Epstein-Barr virus (EBV) is a double-stranded DNA virus and a member of Herpesviridae family. Infection occurs mostly in early childhood and without clinical symptoms. However, infection in teens and young adults may cause infectious mononucleosis. The primary infection is followed by the lifelong latent infection and about 90% of the adult population infected with this virus. EBV has also been linked to a variety of cancers including lymphoma and nasopharyngeal carcinoma. It also can be reactivated in immunocompromised hosts, after transplantation or following therapy in cancer patients

leading to lymphoproliferative disorders. Therefore, viral load monitoring is important in these patients.

#### 4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

#### 5. Kit Contents

The kit contains a manual or quick guide, a CD with Rotor-Gene and StepOne templates and the following reagents:

Label	Content	Quantity
EBV Mix	PCR mix*	360 $\mu$ l
EBV S1	Standard 1: 10,000 copy/ $\mu$ l	150 $\mu$ l
EBV S2	Standard 2: 1,000 copy/ $\mu$ l	150 $\mu$ l
EBV S3	Standard 3: 100 copy/ $\mu$ l	150 $\mu$ l
EBV S4	Standard 4: 10 copy/ $\mu$ l	150 $\mu$ l
Internal Ctrl	Internal control	250 $\mu$ l
Water	PCR Grade Water	200 $\mu$ l

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

#### 6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

## 7. Storage and Stability

The kit components should be shipped and stored at  $-20^{\circ}\text{C}$  and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

## 8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for IVD (in vitro diagnostics) applications.

## 9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and computer accessory
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipetters and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves

## 10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working.**
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

## 11. Specimen, Storage and Transport

Peripheral blood plasma, CSF or tissue should be collected in sterile condition in proper and sterile containers. Samples can be stored at +4°C for a few days or stored at -20°C for up to a few weeks.

## 12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

### 13. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the EBV RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the EBV Mix. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

**Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.**

If the IC is added to EBV Mix, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the PCR Mix before it is added to the tubes. In a successful DNA extraction and PCR test, the IC should generate a CT of 27-34 in the Yellow/VIC Channel.

### 14. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)

- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

## 15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, four for standards and one for the negative control.

**If IC is introduced during the extraction process, pipette 15µl of EBV Mix to each PCR tube.**

**If IC is added to the EBV Mix, add 15ul of the prepared mix (as described in section 13) to each PCR tube.**

**Then add 10ul of extracted DNA, standard, or water to each tube.**

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

*Note: If using Rotor-Gene attach the locking ring.*

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

## 16. Devices and software

EBV RQ kit is designed to work with Rotor-Gene and StepOne.

## 17. Programming Rotor-Gene

*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the CD provided in the kit and double-click on EBV 0.2 or EBV 0.1 according to the tubes used. Program starts. Click on the Start button (Green button on the top menu). On the pop-up

window click Start again and save the program on the desired location.

## 18. Programming StepOne

Open the StepOne software (V 2.\*). On the Set-Up menu, click Template and select the file on the CD provided with the kit. Click on Plate Setup. One negative control, four standards, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment to the desired location. The instrument will start shortly.

## 19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	<b>95°C x 10 min</b>	1
2	<b>95°C x 15 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The EBV Mix contains ROX. Final concentration of ROX in reaction is 300nM.

## 20. Data Analysis: Rotor-Gene

Before analyzing results, make sure that in the sample menu all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be

defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.


Analyze data according to Rotor-Gene manual. Perform quantitative analysis for the **EBV (Green channel)** and qualitative analysis for the **Internal Control (Yellow channel)**. Briefly, click on the analysis menu and then under the Quantitation tab double-click on cycling A. Green.

In the pop-up for Automatic Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK. You may also, simply set the threshold at 0.1. Repeat the above for Cycling A. Yellow, cancel the Automatic Threshold and manually put threshold on 0.1. Figures 1 and 2 represent typical graphs for the Rotor-Gene machine.


**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

Consider following points when analyzing:

- A sample is **Positive** if it is positive in the Green channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in the Green channel while it is positive in the Yellow channel with a sigmoid graph and CT of 27-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of the Green and Yellow channels.

**Fig 1.** Typical EBV graph in Green channel for Rotor-Gene

**Fig 2.** Typical EBV graph in Yellow channel for Rotor-Gene

## 21. Data Analysis: StepOne


Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for the EBV/FAM at 0.1 and 0.05 for the IC/VIC.

Figures 3 and 4 represent typical graphs for the StepOne machine.


**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.**

Consider following points when analyzing:

- A sample is **Positive** if it is positive in the FAM/EBV channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Negative** if it is negative in the FAM/EBV channel while it is positive in the VIC/IC channel with a sigmoid graph and CT of 27-34.
- Results are **inconclusive** and the test should be repeated if a sample is negative in both the FAM/EBV and VIC/IC channels.

**Fig 3.** Typical EBV graph in FAM channel for StepOne

**Fig 4.** Typical EBV graph in VIC channel for StepOne

## 22. Quantitation

The kit provides four quantitation standards with defined titers to generate a standard curve for the quantification of samples, viral load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently, using all four standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ul. To convert the result to copy/ml the following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result(copy/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.

## 23. Linear Range

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 1,000,000 copies/ $\mu\text{l}$  to 1 copy/ $\mu\text{l}$ .

## 24. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned target and showed a limit of detection equal to 0.5 copy/ $\mu\text{l}$ .

## 25. Disposal Method

The contents of the kit do not require any special treatment before disposal and can be directly discarded. Infectious specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

## 26. Technical Support

For technical support, contact us via  
Phone: +98 993-6223241  
email: [info@novingene.com](mailto:info@novingene.com)

## 27. Contact Information

### NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21 88837393

+98 990 11813124




Website: [www.novingene.com](http://www.novingene.com)

## 28. References

- Münz, C. ed., 2015. Epstein Barr virus volume 2: one herpes virus: many diseases (Vol. 391). Springer.
- Kimura, H. and Kwong, Y.L., 2019. EBV viral loads in diagnosis, monitoring, and response assessment. *Frontiers in oncology*, 9, p.62.
- Yin, H., Qu, J., Peng, Q. and Gan, R., 2019. Molecular mechanisms of EBV-driven cell cycle progression and oncogenesis. *Medical microbiology and immunology*, 208, pp.573-583.

- Mackay, Ian M. "Real-time PCR in the microbiology laboratory." *Clinical microbiology and infection* 10, no. 3 (2004): 190-212.

## 29. Symbols

 <b>RUO</b> Research use only	 <b>Manufacturer</b>	 <b>Consult instructions for use</b>
 <b>LOT</b> Lot number	 <b>Content sufficient for &lt;math&gt;\langle \eta \rangle&lt;/math&gt; tests</b>	 <b>Use-by date</b>
 <b>REF</b> Catalogue number	 <b>SN</b> Serial number	 <b>Temperature limit</b>

**For more information and resources please visit our website; [www.novingene.com](http://www.novingene.com)**