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HTLV RQ Kit Manual

For Real-Time PCR Detection and
Quantitation of HTLV-1 provirus

For use with Rotor-Gene or StepOne
Research use only

NG-WI-ASL-13-200
Version 2.0
April 2018


Table of Contents:

1. Introduction.....	2
2. Kit Contents.....	3
3. Storage and Stability.....	3
4. General Precautions.....	3
5. Additionally Required Materials	4
6. Specimen, Storage and Transport.....	4
7. Interfering Substances	5
8. DNA Isolation	5
9. PCR Protocol	5
10. Programming Rotor-Gene	6
11. Programming of StepOne.....	6
12. Programming Other Machines	7
13. Data Analysis: Rotor-Gene.....	7
14. Data Analysis: StepOne.....	9
15. Quantitation	10
16. Linear Range.....	11
17. Sensitivity.....	11

HTLV-I RQ kit is intended for the quantitative detection of HTLV-I provirus DNA with Rotor-Gene or StepOne machines. This kit is for research use only.

1. Introduction

Human T-cell Lymphotropic Virus I (HTLV-I) is the first human retrovirus discovered (1979). Two years later HTLV-II was also discovered. About 15 to 20 million people are infected with these viruses worldwide. HTLV-I is endemic in North of Iran. Transmission is mostly through sexual contact, breastfeeding, Transfusion, Transplant and Intravenous drug use. Within the body, infection spreads slowly through virological synapses from cell to cell and takes 20-50 years to develop symptoms. Infection with HTLV-I may lead to leukemia/lymphoma or myelopathy while HTLV-II causes mild neurologic disorders or chronic pulmonary infection. Disease progress and pathogenesis depend partly on proviral load and monitoring that could be beneficial for patient management.

HTLV-I RQ kit provides a ready-to-use Real-Time PCR system for detection and quantitation of HTLV-I provirus and cellular titer (based on Albumin gene) with Rotor-Gene or StepOne machines. The Albumin gene also monitors extraction quality and prevents false negative results.

Currently Real-Time PCR provides the highest sensitivity and widest dynamic range among other methods for detection of HTLV-I and viral load determination. In this method application of fluorescent probes allows detection of amplified product. Analysis of fluorescent kinetics also leads to quantification of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene and StepOne templates and following reagents:

Label	Content	Quantity
HTLV Mix	PCR Master mix*	480 µl
Albumin Mix	PCR Master mix*	480 µl
HTLV-Alb S1	Standard 1: 100,000 copy/µl	250 µl
HTLV-Alb S2	Standard 2: 10,000 copy/µl	250 µl
HTLV-Alb S3	Standard 3: 1,000 copy/µl	250 µl
HTLV-Alb S4	Standard 4: 100 copy/µl	250 µl
HTLV-Alb S5	Standard 5: 10 copy/µl	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

4. General Precautions

In order to prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction;, b) Reagent preparation where the master-mix is aliquoted into tubes; and c) Reaction preparation area for addition of extracted DNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cold block instead.

5. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

6. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper tubes. We recommend EDTA or citrate tubes. Whole blood should be shipped at +4°C and aliquoted upon receipt. One aliquot should be examined for CBC (see section 9), while other aliquots will be used for DNA extraction and can be stored at +4°C for few days or at -20°C for up to few weeks.

Note: Do not freeze the aliquot specified for CBC!

7. Interfering substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels of

bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

8. DNA Isolation

DNA can be extracted directly from whole blood. Alternatively DNA can be extracted from peripheral blood monocytes (PBMC) collected using ficol.

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

9. PCR Protocol

Each sample should be examined for detection and quantitation of provirus with HTLV-I mix and also with Albumin mix for cell titer. Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider two series of tubes, one for HTLV1 and the other for Albumin. Each of series should include tubes for standards, samples and negative control.

Pipette 20ul of HTLV Mix directly to each tube of series one, followed by adding 5ul of each Standard (HTLV-Alb S1-S5) or sample DNA or control.

Pipette 20ul of Albumin Mix directly to each tube of series two, followed by adding 5ul of each Standard (HTLV-Alb S1-S5) or sample DNA or control.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

10. Programming Rotor-Gene

- Before you start the machine, make sure you have attached the locking ring on the rotor!

On the CD provided with the kit open HTLV Template 0.2 or strip according tubes used. Then, click on Start button (Green button on the top menu). On the pop up window click start again and save program on desired location.

11. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. Positive and Negative controls, standards and few samples are defined. You may change plate set up using right click options (copy, paste, clear). You may also add /remove samples or change sample name on “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment on desired location. Instrument will start shortly.

12. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1 cycle
2	95°C x 15 sec	45 cycles
	60°C x 1 min	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The HTLV and Albumin mixes contain ROX. The final concentration of ROX in reaction is 300nM.

13. Data Analysis: Rotor-Gene

Before analyzing the results, make sure in the sample menu all the standards have been defined as “standard” and the relative concentrations have been entered. **Samples should be named in two separate pages.** Tubes with HTLV mix should be named in “HTLV” page, and tubes with Albumin mix should be named in Albumin page. Define standards in each page relatively. Patient samples should be defined as "unknown" and control as Negative Control" on both pages.

Analyze the data according to manufacturer recommendations. Perform quantitative analysis for **HTLV-I (Green channel)**, and for **Albumin/cell titer (Yellow channel)**. Briefly, click on analysis menu and then, under Quantitation tab double click on “Cycling A. Green”. In the pop up for Automatic Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then, click on OK. Repeat the above for “Cycling A. Yellow”. You may also simply set threshold at 0.1 for both channels.

Refer to figures at the end of the text for typical graphs.

Consider following points when analyzing:

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

- A sample is **positive** if it is positive in both Green and Yellow channels with CTs less than 40. The proviral load can be determined according to section 15.

- A sample is **negative** if it is negative in Green channel (for HTLV-I) while it is positive for Albumin in Yellow channel with a CT within the Albumin standards.
- Results are **inconclusive** if a sample is negative in Yellow channel (for Albumin) and the test should be repeated. Improper DNA extraction or error in test set up could be the cause of such a result.

Figures 1 and 2 represent typical graphs for Rotor-Gene machine.


Figure 1: Typical HTLV-I Graph in Green Channel for Rotor-Gene


Figure 2: Typical Albumin Graph in Yellow Channel for Rotor-Gene

14. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for **HTLV/FAM** at 0.1 and at 0.05 for **Albumin/VIC**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **positive** if it is positive in both FAM and VIC channels with CTs less than 40. The proviral load can be determined according to section 15.
- A sample is **negative** if it is negative in FAM channel (for HTLV-I) while it is positive for Albumin in VIC channel with a CT within the Albumin standards.
- Results are **inconclusive** if a sample is negative in VIC channel (for Albumin) and the test should be repeated.

Improper DNA extraction or error in test set up could be the cause of such a result.

Figures 3 and 4 represent typical graphs for StepOne machine.


Figure 3: Typical HTLV-I Graph in FAM Channel for StepOne


Figure 4: Typical Albumin Graph in VIC Channel for StepOne

15. Quantitation

The kit provides one series of quantitation standards with defined titers for estimation of HTLV-I proviral titer and Albumin gene or cell titer.

If DNA has been extracted from PBMC, the HTLV proviral load is twice the ratio of HTLV/Albumin.

If DNA has been extracted from the whole blood, use the equation below to calculate the proviral load:

$$\text{Proviral Load} = \frac{\text{HTLV titre} \times 2 \times 100\%}{\text{Albumin titer} \times (\text{Lymphocyte \%} + \text{Monocyte \%})}$$

Linear Range

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 1,000,000 copy/ μ l to 10 copy/ μ l.

16. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned target and showed a limit of detection equal to 2 copy/ μ l for HTLV-I and 1 copy/ μ l for Albumin.