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
# **Herpes RQ Kit**

## **Manual**

**For Multiplex Real-Time PCR Detection of CMV,  
HSV-1, HSV-2 , EBV and VZV DNA**

For use with Rotor-Gene  
Research use only

NG-WI-ASL-34-100  
Version 1.1  
Spring 2020

## Table of Contents:

1. Introduction .....	2
2. Kit Contents .....	2
3. Storage and Stability .....	3
4. Additionally Required Materials .....	3
5. General Precautions.....	4
6. Specimen, storage and transport.....	4
7. Interfering substances .....	4
8. Internal control (IC).....	5
9. DNA isolation.....	5
10. PCR Protocol.....	6
11. Programming of the Rotor-Gene .....	6
12. Programming other machines .....	7
13. Data Analysis .....	7
14. Sensitivity.....	11

**Herpes RQ kit** is intended for use with Rotor-Gene and for the detection of Cytomegalovirus (CMV), Herpes Simplex virus-1 (HSV-1), Herpes Simplex virus-2 (HSV-2), Epstein Barr virus (EBV) and Varicella Zoster virus (VZV) DNA. This kit is for research use only.

## 1. Introduction

CMV, HSV-1, HSV-2, EBV and VZV, are all members of Herpesviridae family and have double-stranded DNA genome. These viruses are among common human pathogens and involved in a wide spectrum of diseases including CNS infections.

Herpes RQ kit, provides a fast, sensitive and also cost-effective Multiplex Real-Time PCR system for detection of multiple viruses from Herpesviridae family (CMV, HSV-1, HSV-2, EBV and VZV) in variety of samples. In this method application of fluorescent dye labeled probes allows detection of amplified product in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

This kit is intended for use with Rotor-Gene. This kit also incorporates an Internal Control to identify possible DNA extraction failure or PCR inhibition.

## 2. Kit Contents

The kit contains a manual or quick guide, a CD with Rotor-Gene templates and following reagents:

Label	Content	Quantity
CMV-EBV Mix	PCR mix* for CMV and EBV	360 $\mu$ l
HSV-VZV Mix	PCR mix* for HSV-1, HSV-2 and VZV	360 $\mu$ l
HH Pos Ctrl	Positive control	100 $\mu$ l
HH IC	Internal control	250 $\mu$ l
Water	PCR Grade Water	200 $\mu$ l

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

### 3. Storage and Stability

The kit components should be shipped and stored at  $-20^{\circ}\text{C}$  and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

### 4. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Vortex Mixer
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipments/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

## 5. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components **on ice** completely, mix by flickering followed by a quick spin and **store on crushed ice** after.
- Do not place PCR microtubes on crushed ice. Use cooling blocks instead.

## 6. Specimen, Storage and Transport

Peripheral blood plasma, CSF or tissue should be collected in sterile condition in proper and sterile containers. Samples can be stored at +4°C for few days or stored at -20°C for up to few weeks.

## 7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## 8. Internal Control (IC)

In order to evaluate the possibility of DNA extraction failure and PCR inhibition and prevent false negative results, Herpes RQ kit contains internal control (IC). IC should be used during extraction process and should be added to the mixture of lysis buffer and patient sample. Required volume of internal control is 10% of elution buffer. For instance, if extracted DNA is eluted with 100ul, then 10ul of internal control should be added to the mixture of lysis buffer and patient sample. **Please note that internal control should not be added directly to the patient sample (i.e. before addition of lysis buffer) as it loses its efficiency.**

In a successful DNA extraction and PCR test, IC should generate a CT of 28-32 in Yellow channel with CMV-EBV Mix.

## 9. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAamp UltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

## 10. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin.

Each sample should be evaluated for DNA of CMV, HSV-1, HSV-2, EBV and VZV, in two different series of microtubes. In first series, examining CMV and EBV DNA, consider one tube for each sample plus one for each controls. In second series, examining HSV-1, HSV-2 and VZV DNA, consider one tube for each sample plus one for each controls. Place required number of tubes on cold block.

**In first series of microtubes, pipette 15µl of CMV-EBV Mix to each PCR microtube.**

**In second series of microtubes, pipette 15µl of HSV-VZV Mix to each PCR microtube**

**Then add 10ul of extracted DNA, positive or negative controls to each tube.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

*Note: In Rotor-Gene machine, attach the locking ring.*

## 11. Programming Rotor-Gene

*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the CD provided in the kit and double click on Herpes 0.2 or Herpes 0.1 according to the tubes used. Program starts. Click on Start button (Green button on the top menu). On the pop up window click start again and save program on desired location.

## 12. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	<b>95°C x 10 min</b>	1
2	<b>95°C x 20 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for Green/FAM, Yellow/VIC and Orange/ROX dyes.

## 13. Data Analysis

Analyze data according to manufacturer recommendations. Briefly, click on analysis menu and then under Quantitation tab double click on cycling A. Green. Set the threshold at 0.1. Repeat the above for Yellow and Orange channels and manually put threshold on 0.1. Figures 1, 2 and 3 represent typical graphs for Rotor-Gene.


Fig 1. Typical Controls graph in Green channel for Rotor-Gene



Fig 2. Typical Controls graph in Yellow channel for Rotor-Gene

**Fig 3.** Typical Controls graph in Orange channel for Rotor-Gene

To interpret the results, please note that:

- In **CMV-EBV Mix**, a signal in **Green channel** is due to **EBV**, in **Yellow channel** due to **IC** and in **Orange channel** due to **CMV**.
- In **HSV-VZV Mix**, a signal in **Green channel** is due to **HSV-1**, in **Yellow channel** due to **HSV-2** and in **Orange channel** due to **VZV**.

Table 1 summarizes the above explanation.

Mix	Green Channel	Yellow Channel	Orange Channel
<b>CMV-EBV</b>	<b>EBV</b>	<b>IC</b>	<b>CMV</b>
<b>HSV-VZV</b>	<b>HSV-1</b>	<b>HSV-2</b>	<b>VZV</b>

Table 1: Association of obtained results with evaluated pathogens in two mixes and Green, Yellow and Orange channels of Rotor-Gene

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.**

Consider following points when analyzing:

- A sample is **Positive for EBV**, if it is positive in Green channel with CMV-EBV Mix and a CT of less than 40.
- A sample is **Positive for CMV**, if it is positive in Orange channel with CMV-EBV Mix and a CT of less than 40.
- A sample is **Positive for HSV-1** if it is positive in Green channel with HSV-VZV Mix and a CT of less than 40.
- A sample is **Positive for HSV-2** if it is positive in Yellow channel with HSV-VZV Mix and a CT of less than 40.
- A sample is **Positive for VZV** if it is positive in Orange channel with HSV-VZV Mix and a CT of less than 40.
- A sample is **Negative for CMV and EBV**, if it is negative in Green and Orange channels while it is positive in Yellow channel with a sigmoid graph and a CT of 28-32.
- A sample is **Negative for HSV-1, HSV-2 and VZV**, if it is negative in Green, Yellow and Orange channels while it is positive in Yellow channel of CMV-EBV Mix with a sigmoid graph and CT of 28-32.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in all three channels of Green, Yellow and Orange and the Yellow channel of CMV-EBV Mix.

Interpretation of results is summarized in table 2.

Mix	Green	Yellow	Orange	Result
CMV-EBV	+	+/-	-	Pos: EBV
CMV-EBV	-	+/-	+	Pos: CMV
HSV-VZV	+	-	-	Pos: HSV-1
HSV-VZV	-	+	-	Pos: HSV-2
HSV-VZV	-	-	+	Pos: VZV

Mix	Green	Yellow	Orange	Result
CMV-EBV	-	+	-	Negative
HSV-VZV	-	-	-	

Mix	Green	Yellow	Orange	Result
CMV-EBV	-	-	-	Invalid
HSV-VZV	-	-	-	

Table 2. Interpretation of results in Green, Yellow and Orange channels of Rotor-Gene

## 14. Sensitivity

The analytical detection limit of the kit for each target was assessed with dilution series of the cloned DNA and showed a sensitivity of 0.5 copy/ $\mu$ l for CMV, HSV-1 and HSV-2; and 1 copy/ $\mu$ l for EBV and VZV.