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$$Cut-off - 3.3 < \text{Inconclusive } CT \text{ range} < Cut-off$$

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
# **JAK2 MQ**

# **Kit Manual**

**For Real-Time PCR Detection and  
Quantitation of JAK2 V617F Mutation**

For use with Rotor-Gene or StepOne  
Research use only

NG-WI-ASL-48-100  
Version 1.0  
Winter 2022

## Table of Contents:

1. Introduction .....	2
2. Kit Contents.....	2
3. Storage and Stability .....	3
4. General Precautions.....	3
5. Additionally Required Materials .....	4
6. Specimen, storage and transport.....	4
7. DNA isolation.....	5
8. PCR Protocol.....	5
9. Programming of RotorGene .....	6
10. Programming of StepOne.....	6
11. Programming of other machines .....	7
12. Data Analysis: RotorGene and StepOne .....	7
13. Analytical Sensitivity.....	12

**JAK2 MQ** kit is intended for detection and quantitation of JAK2 V617F mutation in human genomic DNA with Rotor-Gene 6000/Q or StepOne instrument. This kit is for research use only.

## 1. Introduction

JAK2 (Janus Kinase 2) is a Tyrosine Kinase located in cytoplasm with essential role in signaling pathways for cytokines and growth factors. The acquired mutation G1849T replaces valine with phenylalanine (V617F). This substitution results in constitutively active JAK2 which leads to uncontrolled cell proliferation in the absence of growth factors. This mutation is found in the majority of BCR-ABL-negative myeloproliferative disorders (MPDs) and has become a main diagnostic test for polycythemia vera (PV), essential thrombocythemia (ET) and primary myelofibrosis (PMF). JAK2 MQ kit provides a ready-to-use Real-Time PCR system for detection and quantitation of JAK2 V617F mutation with Rotor-Gene or StepOne machines. In this method application of fluorescent dye labeled probes allows detection of amplified product without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

## 2. Kit Contents

The kit contains a manual, quick guide and a CD and following reagents:

Label	Content	Quantity
MJ Mix*	Master mix for JAK2 (V617F)	480 $\mu$ l
WJ Mix*	Master mix for Wild Type	480 $\mu$ l
MJ1	Standard 1 V617F: 100,000 copy/ $\mu$ l	150 $\mu$ l

MJ2	Standard 2 V617F: 10,000 copy/μl	150 μl
MJ3	Standard 3 V617F: 1,000 copy/μl	150 μl
MJ4	Standard 4 V617F: 100 copy/μl	150 μl
MJ5	Standard 5 V617F: 10 copy/μl	150 μl
WJ1	Standard 1 Wild Type: 100,000 copy/μl	150 μl
WJ2	Standard 2 Wild Type: 10,000 copy/μl	150 μl
WJ3	Standard 3 Wild Type: 1,000 copy/μl	150 μl
WJ4	Standard 4 Wild Type: 100 copy/μl	150 μl
WJ5	Standard 5 Wild Type: 10 copy/μl	150 μl
Pos Ctrl 2.5%	Positive Control 2.5%	50 μl
Neg Ctrl	Negative Control	50 μl
Water	PCR Grade Water	200 μl

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

### 3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws especially for PCR Mix more than few times to prevent reduced sensitivity.

### 4. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.

- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and store on **crushed ice while working**.
- Do not place 0.2 ml PCR tubes on crushed ice. Use cold blocks instead.

## 5. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Refrigerated microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 6. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 72 hrs).

DNA can directly be extracted from whole blood. Alternately to increase sensitivity, buffy coat can be used.

Whole blood or buffy coat can be stored at +4°C for three days. Otherwise, should be aliquoted and stored at -70°C which is stable for few months.

## 7. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

To reach the sensitivity of 0.1%, extracted DNA should have a concentration of 10-50 ng/ $\mu$ l. Samples with lower DNA content, result in reduced sensitivity.

## 8. PCR Protocol

Thaw the reagents on ice completely, followed by brief mixing and a quick spin.

Each sample should be examined with both MJ Mix (for V617F) and with WJ Mix (for Wild type allele) in two separate set of reactions. In each set, consider 1 tube for each sample as well as 5 tubes for the 5 standards (MJ or WJ) and 3 tubes for Positive Control, Negative Control and NTC. Use MJ standards for reactions with MJ mix and WJ standards for reactions with WJ mix.

Place required number of tubes on cold block.

**Pipette 20 $\mu$ l of MJ Mix to the first series of tubes and 20 $\mu$ l of WJ Mix to each tube of the second group. Continue by adding 5 $\mu$ l of sample DNA, standards, Positive or Negative Control to each tube.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using RotorGene attach the locking ring.*

## **9. Programming RotorGene**

*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the CD provided in the kit and double click on “JAK2 MQ 0.1” or on “JAK2 MQ 0.2” depending to tubes used.

Click on Start button (Green button on the top menu). On the pop up window click start again and save the run file.

Edit samples on both MJ and WJ pages. Remember that tubes containing MJ Mix should only be named in MJ page and tubes containing WJ Mix should only be named in WJ page.

Make sure in the “Type” column, all the standards have been defined as “standard” and the relative concentrations have been entered. Patient samples should be defined as "unknown" and no template control as "NTC", respectively.

## **10. Programming of StepOne**

Open the StepOne software (V 2.\*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup.

Two targets of MJ and WJ are define in FAM channel. For MJ reactions, 5 standards (MJ1-MJ5) along with Positive and Negative Control, NTC and few samples have been defined. Also, for WJ reactions, 5 standards (WJ1-WJ5) along with Positive and Negative Control, NTC and few samples have been defined. You

may change plate set up using right click options (copy, past, clear). You may also add or remove samples on “Define Targets and Samples” menu. When finished click on Start Run and save the experiment. Instrument will start shortly.

## 11. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	<b>95°C x 3 min</b>	1
2	<b>95°C x 15 sec</b>	50
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for FAM.


Both of MJ Mix and WJ Mix contain ROX with final concentration of 300nM in reaction.

## 12. Data Analysis: RotorGene and StepOne

Analyze the data according to manufacturer recommendations. Perform quantitative analysis for both **MJ** and **WJ** in **Green/FAM channel** with threshold set at 0.1.


Figures 1 and 2 represent typical graphs for StepOne machine. Consider following points when analyzing:

**Fig 1.** Typical MJ graph in FAM channel for StepOne

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.**

**Fig 2.** Typical WJ graph in FAM channel for StepOne

Sample analysis with JAK2 MQ kit can be done in two ways, Quantitative and Qualitative methods.

### A) Quantitative analysis (V617F% calculation)

**Note! To consider a sample valid, total copy number of wild and mutant alleles must be equal to 5,000-30,000 copy/ul.**

V617F% value can be calculated for each patient with below formula.

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- A sample is **Positive** if V617F% is equal or higher than 0.1% and calculated titer could be reported.
- A sample is **Negative** if V617F% is equal or lower than 0.01% and calculated number should be ignored.
- A sample is **Inconclusive** if V617F% is between 0.02% - 0.09%. So it could be Positive or Negative and calculated number should be ignored.

JAK2 V617F%	Results
$\geq 0.1\%$	<b>Positive for V617F</b>
$\leq 0.01\%$	<b>Negative for V617F</b>
$0.02\% \leq \text{Sample} \leq 0.09\%$	<b>Inconclusive</b>

Note: Calculate V617F percentage for JAK2 Neg Ctrl and JAK2 Pos Ctrl. Results must be according to the below table.

Controls	JAK2 V617F%
2.5% Pos Ctrl	1.5% - 4%
Neg Ctrl	$\leq 0.02\%$

## B) Qualitative Analysis

For qualitative detection of JAK2 V617F, in addition to sample, Positive and Negative Controls must be assessed with both MJ and WJ Mix.

First, calculate Cut-off as below:

$$\text{Cut-off} = \Delta \text{CT}_{\text{Neg Ctrl}} = \text{CT}_{\text{MJ Mix}} - \text{CT}_{\text{WJ Mix}}$$

Note! If CT of Negative Control is higher than 45 or it remains negative in Green Channel, consider CT of 45 for it.

- Calculate delta CT of MJ and WJ Mix for each patients sample according to the following equation.

$$\Delta CT \text{ Sample} = CT_{MJ \text{ Mix}} - CT_{WJ \text{ Mix}}$$

- Determine Inconclusive CT range by following formula.

$$\text{Cut-off} < \text{Inconclusive CT range} < \text{Cut-off} - 3.3$$

Consider following points when analyzing:

- A sample is **Positive** for V617F, If sample  $\Delta CT$  is less than Inconclusive CT range.
- A sample is **Negative** for V617F, If sample  $\Delta CT$  is higher than Cut-off.
- A sample is **Inconclusive**, If sample  $\Delta CT$  falls within Inconclusive CT range.

Enter the results in the table below.

Sample	CT <sub>MJ</sub>	CT <sub>WJ</sub>	$\Delta CT$	Inconclusive CT	Evaluation
Pos Ctrl					
Neg Ctrl					
Sample					
Sample					
Sample					

As an example, If CT of Neg Ctrl with WJ is 24.8 in Green/FAM channel and 37.2 with MJ Mix, the  $\Delta CT$  or Cut-off is 12.4. and Inconclusive CT range is  $(12.4 - 3.3 = 9.1)$  9.1 to 12.4. results are summarized in below table.

Sample	CT <sub>MJ</sub>	CT <sub>WJ</sub>	ΔCT	Inconclusive CT	Evaluation
Neg Ctrl	37.2	24.8	12.4	9.1-12.4	Negative
Pos Ctrl	29.2	23.4	5.8	5.8 < 9.1	Positive
Sample	41.5	28.5	13	13 > 9.1	Negative
Sample	30.7	29.0	1.7	1.7 < 9.1	Positive
Sample	39.5	28.7	10.8	9.1 < 10.8 < 12.4	Inconclusive

### 13. Analytical Sensitivity

Analytical sensitivity of the kit has been estimated about 0.1% for JAK2 V617F. To achieve this sensitivity, sample must include 5,000-30,000 copy/μl mutant and wild type alleles.