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ایمیل: info@novingene.com

وبسایت: www.novingene.ir

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- Jones, D. (2010). Neoplastic Hematopathology. Springer Science & Business Media.
- Yoo, J.H., Park, T.S., Maeng, H.Y., Sun, Y.K., Kim, Y.A., Kie, J.H., Cho, E.H., Song, J., Lee, K.A., Suh, B. and Choi, J.R., 2009. JAK2 V617F/C618R mutation in a patient with polycythemia vera: a case study and review of the literature. Cancer genetics and cytogenetics, 189(1), pp.43-47.
- Levine, R.L., Pardanani, A., Tefferi, A. and Gilliland, D.G., 2007. Role of JAK2 in the pathogenesis and therapy of

myeloproliferative disorders. Nature reviews cancer, 7(9), pp.673-683.

- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.

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JAK2 RQ Kit Manual

Spring 2024, Version 5.0

For Real-Time PCR Detection of JAK2 V617F Mutation
For use with Rotor-Gene or StepOne
For Research use only

 24 (Cat# JAK2RQ24)

 48 (Cat# JAK2RQ48)

 96 (Cat# JAK2RQ96)

 NG-WI-ASL-10-500

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use.....	3
3. Background Information	3
4. Test Principle	3
5. Kit Contents	4
6. Packaging models.....	4
7. Storage and Stability	4
8. Product Use Limitations	4
9. Additionally Required Materials.....	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Interfering Substances	6
13. DNA Isolation	6
14. Internal Control (IC).....	7
15. PCR Protocol	7
16. Devices and software	7
17. Programming Rotor-Gene.....	7
18. Programming StepOne	8
19. Programming Other Machines	8

20. Data Analysis: Rotor-Gene.....	8
21. Data Analysis: StepOne	11
22. Sensitivity.....	14
23. Disposal Method	14
24. Technical Support	14
25. Contact Information.....	14
26. References.....	15
27. Symbols	15

1. Introduction

JAK2 RQ kit provides a ready-to-use Real-Time Time polymerase chain reaction (PCR) test designed for detecting JAK2 V617F mutation. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probes for detecting a housekeeping gene sequence as an Internal Control (IC). IC prevents false negative results due to extraction or presence of PCR inhibitors.

This kit is intended for Research Use Only!

2. Intended Use

JAK2 RQ kit is intended for detecting JAK2 V617F mutation in DNA samples extracted from human. Detection is achieved using Real-Time PCR and is compatible with Rotor-Gene or StepOne machines.

3. Background Information

JAK2 (Janus Kinase 2) is a Tyrosine Kinase located in cytoplasm with essential role in signaling pathways for cytokines and growth factors. The acquired mutation G1849T replaces valine with phenylalanine (V617F). This substitution results in constitutively active JAK2 which leads to uncontrolled cell proliferation in the absence of growth factors. This mutation is found in the majority of BCR-ABL-negative myeloproliferative disorders (MPDs) and has become a main diagnostic test for polycythemia vera (PV), essential thrombocythemia (ET) and primary myelofibrosis (PMF).

4. Test Principle

The target sequence is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through

JAK2 RQ (V5.0)

fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual, a CD with Rotor-Gene and StepOne templates and the following reagents:

Label	Content	Quantity
JAK2 Mix	PCR Master mix*	480 µl
Pos Ctrl 5%	Positive Control 5%	100 µl
Pos Ctrl 0.1%	Positive Control 0.1%	100 µl
Negative Ctrl	Negative Control	100 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The User manual should be strictly followed, and any modification will invalidate the results.

- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for in vitro diagnostics.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and computer accessory
- Tabletop microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.5ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

To prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the Master Mix is aliquoted into tubes, and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.

- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

11. Specimen, Storage and Transport

Peripheral blood is the suitable sample for JAK2 test. We recommend EDTA or citrate as anticoagulant.

DNA can directly be extracted from whole blood. Alternately to increase sensitivity buffy coat can be used.

Whole blood or buffy coat can be shipped and stored at +4°C for a few days or can be aliquoted and stored at -20°C which is stable for a few months.

12. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes and samples of heparinized patients must not be used.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

13. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany).

To reach the sensitivity of 0.1%, extracted DNA should have a concentration of 10-50 ng/ul.

14. Internal Control (IC)

To examine DNA extraction quality as well as the presence of PCR inhibitors and to prevent false negative results, primers and probe for an *Internal Control* (a housekeeping gene) is included in the PCR Master Mix. Internal control should generate a CT of 22-25 in Rotor-Gene and a CT of 24-28 in StepOne. See analysis section for more details.

15. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, plus three for Positive control, Negative control and NTC.

Pipette 20ul of JAK2 Mix directly, to each tube followed by adding 5ul of control, or isolated DNA.

Cap the tubes and visually inspect to make ensure all are capped securely. Place the tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

16. Devices and software

JAK2 RQ kit is designed to work with Rotor-Gene and StepOne.

17. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the CD provided in the kit and double-click on "JAK2 0.2" or on "JAK2 strip" depending on the tubes used. Program starts. Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

18. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu, click Template and select the file on the CD provided with the kit. Click on Plate Setup. One negative control, four standards, and a few samples are defined. You may change plate setup using right-click options (copy, paste, clear). You may also add /remove samples or change the sample name on the "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment to the desired location. The instrument will start shortly.

19. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 15 sec	50
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC/HEX dyes. JAK2 Mix contains 300nM ROX in final reaction.

20. Data Analysis: Rotor-Gene

Analyze data according to manufacturer recommendations. Perform qualitative analysis for **JAK2 (Green channel)** and **Internal Control (Yellow channel)**. Briefly, click on "Analysis" menu and then under "Quantitation" tab, double click on "Cycling A. Green". Close the pop up for Automatic Threshold and set the threshold on 0.1 and above background or Negative control. Repeat the above for "Cycling A. Yellow" and manually set threshold on 0.1. Refer to figures 1 and 2 for typical graphs.

JAK2 RQ (v5.0)


Fig 1. Typical JAK2 graph in Green channel for Rotor-Gene


Fig 2. Typical JAK2 graph in Yellow channel for Rotor-Gene

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

For detection of JAK2 V617F, in addition to sample, 0.1% Positive Control and Negative Control must be assessed too. Calculate ΔCT for each sample including 0.1% Pos Ctrl and Neg Ctrl and patient sample according to the below formula.

$$\Delta CT = CT_{\text{Green Ch}} - CT_{\text{Yellow Ch}}$$

Note! If CT of Negative Control is higher than 45 or it remains negative in the Green Channel, use CT of 45 for calculation.

Consider the following points for analyzing:

- A sample is **Positive** for V617F, when sample ΔCT is less than 0.1% Pos Ctrl ΔCT .
- A sample is **Negative** for V617F, when sample ΔCT is higher than Neg Ctrl ΔCT .
- A sample is **Inconclusive**, when sample ΔCT is higher than 0.1% Pos Ctrl ΔCT and less than Neg Ctrl ΔCT .

Above points are summarized in the following Table.

ΔCT	Result
$\Delta CT_{\text{Sample}} < \Delta CT_{\text{0.1\% Pos Ctrl}}$	Positive for V617F
$\Delta CT_{\text{Sample}} > \Delta CT_{\text{Neg Ctrl}}$	Negative for V617F
$\Delta CT_{\text{0.1\% Pos Ctrl}} < \Delta CT_{\text{Sample}} < \Delta CT_{\text{Neg Ctrl}}$	Inconclusive

As an example, If CT of Neg Ctrl is 47 in Green channel and 24 in the Yellow channel, the ΔCT is 23. And, If CT of 0.1% Pos Ctrl is 37 in the Green channel and 24 in the Yellow channel, the ΔCT is 13. For patient sample with the CT of 36 in the Green channel and 24 in the Yellow channel, ΔCT is 12 and less than 0.1% Pos Ctrl ΔCT (13) so, patient sample is considered Positive for V617F. The example is summarized in below Table.

Sample	Green CT	Yellow CT	ΔCT
Neg Ctrl	47	24	47-24=23
0.1% Pos Ctrl	37	24	37-24=13
Patient sample	36	24	36-24=12

21. Data Analysis: StepOne

Analyze the data according to the manufacturer's recommendations. Briefly, click on the "Analyze" and set the threshold at 0.1 for **JAK2 (the FAM channel)** and 0.05 for **Internal Control (the VIC channel)**.

Refer to figures 3 and 4 for typical graphs.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.


Fig 3. Typical JAK2 graph in FAM channel for StepOne


Fig 4. Typical JAK2 graph in VIC channel for StepOne

For detection of JAK2 V617F, in addition to sample, 0.1% Positive Control and Negative Control must be assessed too. Calculate ΔCT for each sample including 0.1% Pos Ctrl and Neg Ctrl and patient sample according to the below formula.

$$\Delta CT = CT_{FAM Ch} - CT_{VIC Ch}$$

Note! If CT of Negative Control is higher than 45 or it remains negative in Green Channel, use CT of 45 for calculation.

Consider following points for analyzing:

- A sample is **Positive** for V617F, when sample ΔCT is less than 0.1% Pos Ctrl ΔCT .
- A sample is **Negative** for V617F, when sample ΔCT is higher than Neg Ctrl ΔCT .
- A sample is **Inconclusive**, when sample ΔCT is higher than 0.1% Pos Ctrl ΔCT and less than Neg Ctrl ΔCT .

Above points are summarized in the following table.

ΔCT	Result
$\Delta CT \text{ Sample} < \Delta CT \text{ 0.1\% Pos Ctrl}$	Positive for V617F
$\Delta CT \text{ Sample} > \Delta CT \text{ Neg Ctrl}$	Negative for V617F
$\Delta CT \text{ 0.1\% Pos Ctrl} < \Delta CT \text{ Sample} < \Delta CT \text{ Neg Ctrl}$	Inconclusive

As an example, If CT of Neg Ctrl is 47 in FAM channel and 28 in VIC channel, the ΔCT is 19. And, If CT of 0.1% Pos Ctrl is 39 in FAM channel and 28 in VIC channel, the ΔCT is 11. For patient sample with the CT of 39 in FAM channel and 28 in VIC channel, ΔCT is 8 and less than 0.1% Pos Ctrl ΔCT (11) so, patient sample is considered Positive for V617F. The example is summarized in below Table.

Sample	FAM CT	VIC CT	Δ CT
Neg Ctrl	47	28	47-28=19
0.1% Pos Ctrl	39	28	39-28=11
Patient sample	35	27	35-27=8

22. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the positive target in negative DNA and showed a limit of detection equal to 0.1% or 1 in 1,000. In other words, if in 1,000 cells only one is positive, it will be detected in 95% cases. To meet this level of sensitivity, extracted sample should contain 10-50 ng/ μ l DNA.

23. Disposal Method

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

24. Technical Support

For technical support, contact us via phone at +98 993-6223241 and email: info@novingene.com

25. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124



Email: info@novingene.com

Website: www.novingene.com

26. References

- Jones, D. (2010). Neoplastic Hematopathology. Springer Science & Business Media.
- Yoo, J.H., Park, T.S., Maeng, H.Y., Sun, Y.K., Kim, Y.A., Kie, J.H., Cho, E.H., Song, J., Lee, K.A., Suh, B. and Choi, J.R., 2009. JAK2 V617F/C618R mutation in a patient with polycythemia vera: a case study and review of the literature. Cancer genetics and cytogenetics, 189(1), pp.43-47.
- Levine, R.L., Pardananani, A., Tefferi, A. and Gilliland, D.G., 2007. Role of JAK2 in the pathogenesis and therapy of myeloproliferative disorders. Nature reviews cancer, 7(9), pp.673-683.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.

27. Symbols

 RUO Research use only	 Manufacturer	 Consult instructions for use
 LOT Lot number	 Content sufficient for <n> tests	 Use-by date
 REF Catalogue number	 SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com