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STI-3 RQ Kit Manual

Summer 2024, Version 1.0

For Real-Time PCR Detection of *Chlamydia trachomatis*,
Neisseria gonorrhoeae, and *Mycoplasma genitalium*
For use with Rotor-Gene
For Research use only

 24 (Cat#STI3RQ24)

 48 (Cat#STI3RQ48)

 96 (Cat#STI3RQ96)

 NG-WI-ASL-62-100

RUO

NovinGene ParsVira

No. 5, Iranshahr St, Tehran, Iran 1581633336.


Table of Contents

1. Introduction	3
2. Intended Use	3
3. Background Information	3
4. Test Principle	4
5. Kit Contents	4
6. Packaging models	4
7. Storage and Stability	4
8. Product Use Limitations	5
9. Additionally Required Materials.....	5
10. General Precautions	5
11. Specimen, Storage and Transport	6
12. Internal Control (IC)	6
13. DNA Isolation	7
14. PCR Protocol	7
15. Devices and software.....	8
16. Programming of the Rotor-Gene.....	8
17. Programming Other Machines	8
18. Data Analysis: Rotor-Gene	8
19. Sensitivity.....	12

20. Disposal Method	12
21. Technical Support.....	13
22. Contact Information.....	13
23. References	13
24. Symbols.....	14

1. Introduction

STI-3 RQ kit provides a ready-to-use Real-Time polymerase chain reaction (PCR) test designed for detecting DNA of the following bacteria: *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, and *Mycoplasma genitalium*. All required reagents are included in the PCR Mix provided in the kit. This Mix also contains a different series of primers and probe for detecting a synthetic DNA sequence. The kit supplies this synthetic DNA sequence as an Internal Control (IC). The IC can be used either during DNA extraction or in the PCR reaction to prevent false negative results due to failure in the above steps.

This kit is intended for Research Use Only!

2. Intended Use

STI-3 RQ kit is intended for detecting DNA of following bacteria, *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, and *Mycoplasma genitalium*. Detection is achieved using Multiplex Real-Time PCR system with Rotor-Gene machine.

3. Background Information

Sexually Transmitted Infections (STIs) represent a significant global public health concern. They can be caused by bacteria, viruses or parasites. Most common bacteria transmitted through sexual contact including; *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, and *Mycoplasma genitalium*. Infected persons may exhibit various symptoms, although in many cases infected person may remain asymptomatic. Delayed diagnosis and treatment may lead to serious and sometimes permanent health complications such as pelvic inflammatory diseases and infertility.

4. Test Principle

The pathogen is detected using PCR, where primers specific to the target genome amplify its unique sequence. Real-Time PCR facilitates the detection of the amplified product through fluorescent-labeled probes. Therefore, monitoring fluorescence provides a means for detecting the target without requiring post-amplification analysis. This eliminates the possibility of PCR product contamination.

5. Kit Contents

The kit contains a manual a CD with Rotor-Gene templates and the following reagents:

Label	Content	Quantity
STI Mix A	PCR Master mix* for Chlamydia trachomatis, Neisseria gonorrhoeae, and Mycoplasma genitalium	360 µl
STI-3 Pos Ctrl	Positive Control	150 µl
Internal Ctrl	Internal Control*	250 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

6. Packaging models

The kit is available in 24, 48, and 96 reactions of 25 microliters.

7. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiration date mentioned. Avoid repeated freeze-thaw more than three times to prevent reduced sensitivity.

8. Product Use Limitations

- This kit is intended to be used only by specially instructed and trained personnel.
- The user manual should be strictly followed, and any modification will invalidate the results.
- The kit and its contents should not be used past the expiration date on the package.
- The kit and its contents should not be used if there is any sign of pink or red color on the Warm Mark label.
- This kit is for Research use only and is not validated for in vitro diagnostics.

9. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and computer accessory
- Tabletop microtube centrifuge
- Vortex Mixer
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipments/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

10. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the PCR mix is aliquoted into tubes, and

- c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
 - Thaw on ice kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.
 - Keep PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.
 - Do not place 0.2ml PCR tubes on crushed ice. Use cold block instead.

11. Specimen, Storage and Transport

Proper samples to test is urine (20ml of voided urine), cytobrush, swabs include: vaginal, cervical and male urethral. Samples can be stored at $2-8^{\circ}\text{C}$ for a few hours or at -20°C or lower for a few months.

12. Internal Control (IC)

To assess the possibility of DNA extraction failure and PCR inhibition, and to prevent false negative results, the STI-3 RQ kit contains an internal control (IC). This IC can be used during the extraction process or added directly to the STI Mix A. To monitor both DNA extraction and PCR reaction, the IC should be added to the mixture of lysis buffer and patient sample during extraction. The required volume of IC is 10% of the elution buffer. For instance, if the extracted DNA is eluted with 100ul, then 10ul of IC should be added to the mixture of lysis buffer and patient sample.

Please note that the IC should not be added directly to the patient sample (i.e, before the addition of lysis buffer) as it loses its efficiency.

If the IC is added to STI Mix A, only PCR inhibition can be monitored. For this purpose, 1ul of the IC should be added to each reaction. For example, for 10 PCR reactions, 10ul of the IC should be added to 150ul of the STI Mix A before it is added to the tubes. In a successful DNA extraction and PCR test, the IC should generate a CT of 27-32 in the Yellow/VIC Channel.

13. DNA Isolation

DNA isolation from bacteria can be performed using different kits from various manufacturers. We recommend the following:

- QIAamp DNA Blood Mini Kit (Cat# 51104, Qiagen, Germany)
- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)

14. PCR Protocol

Thaw the reagents on ice completely, followed by a brief mixing and a quick spin. Place the required number of tubes on a cold block. Consider one tube for each sample, plus one for positive sample and one for water.

If IC is introduced during the extraction process, pipette 15µl of the [STI Mix A Mix](#) to each PCR tube.

If the IC is added to the [STI Mix A Mix](#), add 15ul of the [prepared mix](#) (as described in section 12) to each PCR tube.

Then add 10ul of isolated DNA, [Pos Ctrl](#), or water to each tube.

Cap the tubes and visually inspect to ensure all are capped securely. Place the tubes in the machine.

Note: If using Rotor-Gene attach the locking ring.

15. Devices and software

STI-3 RQ kit is designed to work with Rotor-Gene.

16. Programming of the Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the CD provided in the kit and double-click on Respiratory template "STI3 0.1" or "STI3 0.2" depending on the tubes used.

Program starts. Click on the Start button (Green button on the top menu). On the pop-up window click Start again and save the program on the desired location.

17. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table (please note that following program is not applicable to Rotor-Gene):

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM, VIC, ROX and Cy5 dyes.

18. Data Analysis: Rotor-Gene

To analyze data briefly, click on the analysis menu and then under Quantitation tab double click on cycling A. Green. Manually put threshold, at 0.1. Repeat the above for Yellow, Red and Orange Channels.

Figures 1, 2, 3 and 4 represent typical graphs for Rotor-Gene. To interpret the results, please note that:

A Signal in the **Green** channel indicates **Chlamydia trachomatis**, **Orange** channel indicates **Mycoplasma genitalium** and **Red** channel is due to **Neisseria gonorrhoeae**, and a signal in the **Yellow** channel is due to **Internal Control**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of a sigmoid graph and log phase, the sample is considered Negative, and CT, if present, is not reliable.

Consider the following points when analyzing:

- A sample is **Positive** for Chlamydia trachomatis if it is positive in the **Green** channel with sigmoid graphs and CT of less than 40 and also positive in the **Yellow** channel with a CT of 20-30.
- A sample is **Positive** for Mycoplasma genitalium if it is positive in the **Orange** channel with sigmoid graphs and CT of less than 40 and also positive in the **Yellow** channel with a CT of 20-30.
- A sample is **Positive** for Neisseria gonorrhoeae if it is positive in the **Red** channel with sigmoid graphs and CT of less than 40 and also positive in the **Yellow** channel with a CT of 20-30.
- A sample is **Negative** for Chlamydia trachomatis, Mycoplasma Genitalium and Neisseria gonorrhoeae if it is negative in the **Green** and the **Yellow** and the **Red** channels while it is positive in the **Yellow** channel with a sigmoid graph and CT of 20-30.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in **Yellow** channel.

The interpretation of results is summarized in bellow Table1.

Green	Yellow	Red	Orange	Result
+	+	-	-	Pos for Chlamydia trachomatis
-	+	-	+	Pos for Mycoplasma genitalium
-	+	+	-	Pos for Neisseria gonorrhoeae
-	+	-	-	Negative
- / +	-	- / +	- / +	Inconclusive


Fig 1. Typical Controls graph in Green channel for Rotor-Gene (Chlamydia trachomatis)


Fig 2. Typical Controls graph in Yellow channel for Rotor-Gene (Internal Control)


Fig 3. Typical Controls graph in Red channel for Rotor-Gene (*Neisseria gonorrhoeae*)


Fig 4. Typical Controls graph in Orange channel for Rotor-Gene (Mycoplasma genitalium)

19. Sensitivity

The analytical detection limit of the kit was assessed with a cloned targets genomes and showed a limit of detection equal to 2.8 copies/ μ l for *Chlamydia trachomatis* and 3.8 copies/ μ l for *Neisseria gonorrhoeae*, and 4 copies/ μ l for *Mycoplasma genitalium*.

20. Disposal Method

The contents of the kit do not require any special method of disposal and can be directly discarded. But infectious laboratory specimens should be maintained in 5% Sodium Hypochlorite overnight and then discarded.

21. Technical Support

For technical support, contact us via phone at +98 993-6223241 and email: info@novingene.com

22. Contact Information

NovinGene ParsVira

Address: No. 5, Iranshahr St, Tehran, Iran 1581633336.

Tel: +98 21-88837393

+98 990 -1813124

Email: info@novingene.com

Website: www.novingene.com

23. References

- Mahon, C.R. and Lehman, D.C. (2023). Textbook of Diagnostic Microbiology - E-Book. Elsevier Health Sciences.
- CDC (2023). STI Treatment Guidelines. [online] Centers for Disease Control and Prevention. Available at: <https://cdc.gov/std/treatment-guidelines/default.htm> [Accessed 21 May 2024].
- Buder, S., Schöfer, H., Meyer, T., Bremer, V., Kohl, P.K., Skaletz-Rorowski, A. and Brockmeyer, N. (2019). Bacterial sexually transmitted infections. JDDG: Journal der Deutschen Dermatologischen Gesellschaft, 17(3), pp.287–315. doi:<https://doi.org/10.1111/ddg.13804>.
- Mackay IM. Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. 2004; 10 (3): 190 – 212.

24. Symbols

RUO Research use only	 Manufacturer	 Consult instructions for use
LOT Lot number	 Content sufficient for <n> tests	 Use-by date
REF Catalogue number	SN Serial number	 Temperature limit

For more information and resources please visit our website; www.novingene.com