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EGFR RQ Kit Manual

**For Real-Time PCR
Detection of EGFR Mutations**

Research use only

NG-WI-ASL-42-102
Version 1.2
Autumn 2021


Table of Contents:

1. Introduction.....	2
2. Kit Contents.....	4
3. Storage and Stability	4
4. General Precautions.....	5
5. Additionally Required Items.....	5
6. Specimen	5
7. DNA Extraction	6
8. Protocol: DNA Sample Assessment.....	6
9. Programming Real-Time PCR.....	7
10. Analysis: Sample Assessment.....	7
11. Protocol: Detection of EGFR mutation.....	9
12. EGFR Mutation Detection Analysis.....	11
13. Analytical Sensitivity	19

EGFR RQ kit is intended for the detection of 34 somatic mutations in the EGFR gene. This kit is designed for use with Real-Time PCR machine. This kit is for research use only!

1- Introduction

Epidermal Growth Factor Receptor (EGFR), is a tyrosine kinase receptor and is considered as an oncogene. EGFR is involved in regulation of cellular proliferation, differentiation and survival. Mutations in EGFR exons 18, 19, 20 and 21 are associated with the development of different human cancers specially None Small Cell Lung Carcinoma (NSCLC), and glioblastoma. Since, the choice of anti-EGFR therapies are highly dependent on the EGFR mutations, it is essential that patients are tested for these mutations.

This kit provides ready to use reagents for detection of 34 mutations. Detected mutations include three point mutations in exon 18 (G719A, G719S and G719C without differentiating them), 24 deletions in exon 19 (without differentiation between them), three insertion (without differentiation between them) and two point mutation in exon 20 (S768I and T790M) and two point mutations in exon 21 (L858R, L861Q). For more information see table1.

Exon	Mutation	Cosmic* ID	Base change
18	G719A	6239	2156G>C
	G719S	6252	2155G>A
	G719C	6253	2155G>T
19	Deletions	6210	2240_2251 del 12
		6218	2239_2247 del 9

		6220	2238_2255 del 18
		6223	2235_2249 del 15
		6225	2236_2250 del 15
		6254	2239_2253 del 15
		6255	2239_2256 del 18
		12367	2237_2254 del 18
		12369	2240_2254 del 15
		12370	2240_2257 del 18
		12382	2239_2248TTAAGAGAAG>C
		12383	2239_2251>C
		12384	2237_2255>T
		12385	2235-2255delinsAAT
		12387	2239_2258>CA
		12403	2239-2256 delinsCCA
		12419	2238_2252>GCA
		12422	2238_2248>GC
		12678	2237_2251 del 15
		12728	2236_2253 del 18
		13550	2235-2248delinsAATTC
		13551	2235_2252>AAT
		13552	2235-2251delinsAATTC
		23571	Del 2238-2252
20	S768I	6241	2303G>T
	Insertion	12376	2307_2308insGCCAGCGTG
		12378	2310_2311insGGT
		12377	2319_2320insCAC
T790M	6240	2369C>T	
21	L858R	6224	2573T>G
	L861Q	6213	2582T>A

Table1: EGFR mutations detected by NovinGene kit.

* COSMIC:Catalogue of somatic mutations in cancer: <http://cancer.sanger.ac.uk/>.

This kit is based on based on Real-Time PCR technology. In this method application of fluorescent probes and analysis of

fluorescent kinetics allows detection of amplified product without requiring post-amplification analysis. Therefore, test is accomplished much faster and possibility of contamination with the PCR product is reduced sharply. This kit also includes internal control for preventing false negative results due to failure in DNA extraction or PCR reaction.

2- Kit Contents

The kit contains a manual and CD and following reagents:

Lable	Content	Quantity	Volume
EGFR Ctrl Mix	PCR Mix for quality control	2	480 µl
G719X Mix	PCR Mix for G719A/G719SIG719C	1	480 µl
19Del Mix	PCR Mix for exon 19 Deletions	1	480 µl
20Ins Mix	PCR Mix for exon 20 insertions	1	480 µl
S768I Mix	PCR Mix for S768I	1	480 µl
T790M Mix	PCR Mix for T790M	1	480 µl
L858R Mix	PCR Mix for L858R	1	480 µ
L861Q Mix	PCR Mix for L861Q	1	480 µ
EGFR Pos	Positive control	1	250 µl
EGFR Neg	Negative control	1	250 µl
Water	PCR Grade Water	1	200 µl

3- Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws especially for PCR Mix more than few times to prevent reduced sensitivity.

4- General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

5- Additionally Required Items

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Adjustable pipettors and nuclease free filtered tips
- Vortex Mixer
- DNA extraction kit and required equipment/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

6- Specimen

Sections of formalin-fixed tissue (paraffin block) or unfixed tissue can be examined. Note that tumor tissue is heterogeneous with

uneven distribution of mutant cells; therefore different parts of tissue may produce different results.

7- DNA Extraction

DNA extraction can be performed with different kits from various manufacturers. We recommend using:

- QIAamp Fast DNA Tissue Kit (Cat. no. 51404, Qiagen GmbH, Hilden, Germany)
- Roche; High Pure FFPE DNA Isolation Kit (Cat. No. 06 650 767 001, Roche Diagnostics GmbH, Mannheim, Germany)

Note that, the minimum volume of DNA required for the test is 40 microliters of 5-30 ng/ul. We recommend extracting 100 microliters or higher volume of DNA, enough to repeat the test if necessary.

8- Protocol: DNA Sample Assessment

Before examining a sample for EGFR mutations, quality of DNA should be assessed. If the results are within the desired range, then the second test for detecting EGFR mutations will be performed. To qualify DNA extraction, follow steps below.

Note! Check the volume of extracted DNA sample and make sure it is more than 40 microliters and preferably about 100 microliters. Volumes less than 40 microliters are not enough to proceed.

First, thaw the **EGFR Ctrl Mix** on ice completely, mix by inversion followed by a quick spin and store on crushed ice after. Place the required number of microtubes on cold block including one for each sample, plus two for negative control and water.

Pipette 20µl of EGFR Ctrl Mix to each microtube. Continue by adding 5µl of DNA sample, Negative control and Water to each tube.

Cap the tubes and inspect visually to make sure all are capped securely. Place tubes in the machine.

9- Programming Real-time PCR

Rotor-Gene or StepOne could be setup using templates provided in the CD. Templates may also be downloaded from our website www.novingene.com. For Rotor-Gene, select “EGFR 0.2” or “EGFR 0.1” according to the PCR tubes.

You may also set up the machine as below.

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM/Green and VIC/Yellow dyes/channels.

All mixes contain ROX with final concentration of 300nM in the reaction.

10-Analysis: Sample Assessment

Data Analysis is performed in two steps. First, the test validity is assessed and then DNA sample analysis. As the initial step, results of negative control will be evaluated.

Perform quantitative analysis for both **DNA sample (FAM/Green channel)** and **internal control (VIC/Yellow channel)** with the threshold set at 0.05 for both channels.

Note that, a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample can not be considered as positive (even CT present) result is not reliable.

If above are met, then the test results are valid and may proceed to DNA sample assessment.

Test is valid only if:

- 1) Water sample is **Negative** in FAM/Green channel and it is **Positive** in VIC/Yellow channel with CT of 27-33.
- 2) Negative controls is **Positive** in VIC/Yellow channel with CT of 27-33.
- 3) Negative control is **Positive** in FAM/Green channel with CT of 23-28.

Expected results for a valid test are summarized in Table2.

Sample	VIC/Yellow	FAM/Green
NTC	Pos (CT: 27-33)	Neg
Negative Ctrl	Pos (CT: 27-33)	Pos (CT: 23-28)

Table2. Expected results for valid Sample Assessment.

B) DNA sample Quality Analysis:

DNA sample is qualified only if sample is positive in FAM/Green channel with CT of 22-30 and positive in VIC/Yellow

channel with CT of 27-33. This sample can be further examined for EGFR mutations.

Sample should be diluted with water, if a sample is positive in FAM/Green channel with CT of less than 22 and positive in VIC/Yellow channel with CT of 27-33. Dilute the sample to reach CT of 22-27. By 2X dilution of sample, CT is increased one cycle. After diluting the sample, proceed to the EGFR mutation test.

Result is invalid, if a sample in FAM/Green channel has CT above 30 and in VIC/Yellow channel CT of 27-33. Also result is invalid if a sample in VIC/Yellow channel has CT above 33 and in FAM/Green channel has CT of 22-30. In this case, repeating DNA extraction is recommended.

DNA sample Quality Analysis is summarized in Table 3.

	FAM/Green	VIC/Yellow	Result
1	+ CT: 22-30	+ CT: 27-33	Valid
2	+ CT<22	+ CT: 27-33	Sample dilution
3	+ CT>30	+ CT: 27-33	Invalid
4	+ CT: 22-30	+ CT>33	Invalid

Table3. DNA Quality Analysis

11-Protocol: Detection of EGFR mutations

To detect EGFR mutations, each sample should be examined with eight mixes, one mix for DNA quantitation and seven mixes for detection of 34 mutations. Therefore, to examine only one

EGFR RQ (V1.2)

sample, 32 microtubes are required; 8 for the Positive control, 8 for the Negative control and 8 for NTC and 8 for the sample. Respectively for each extra sample, 8 microtubes will be added. Therefore, for two samples, 40 microtubes and for three samples 48 microtubes are required. Figure1 shows tube setup for three samples.

To start, Place required number of tubes on cold block organized in series of six each.

Pipette 20 μ l of EGFR Ctrl Mix to each tube in the first series.

Pipette 20 μ l of G719X Mix to each tube in the second series.

Pipette 20 μ l of 19Del Mix to each tube in the third series.

Pipette 20 μ l of 20Ins Mix to each tube in the fourth series.

Pipette 20 μ l of S768I Mix to each tube in the fifth series.

Pipette 20 μ l of T790M Mix to each tube in the sixth series.

Pipette 20 μ l of L858R Mix to each tube in the sixth series.

Pipette 20 μ l of L861Q Mix to each tube in the sixth series.

Continue by adding 5 μ l of Positive control, Negative Control, water and extracted DNA to each tube. Consider the first and second tube in each series for positive control and Negative control and the third for water/NTC. The next tubes would be for samples.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine and program it according to the section 9.

Note: if using Rotor-Gene attach the locking ring too.


Fig1. Tubes setup for Mixes and samples.

12-EGFR Mutation Detection Analysis

Before analyzing results for EGFR mutations, test validity should be verified. To do so, results of positive and negative controls and samples with Control Mix will be examined. If the test is valid, then may proceed to EGFR mutation analysis.

To analyze results, set threshold at 0.05 for FAM/Green and VIC/Yellow channel.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable. In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

A) Test is valid only if:

- 1) NTC is Negative in FAM/Green channel with all Mixes, and
- 2) Positive Ctrl, Negative Ctrl and NTC are Positive in VIC/Yellow channel with each Mix with a CT of 27-33, and
- 3) Positive control in FAM/Green channel with Control Mix with a CT of 22-27, and
- 4) Positive control is Positive in FAM/Green channel with all mixes with of 27-35,
- 5) Sample is Positive in VIC/Yellow channel with Control Mix with a CT of 27-33,
- 6) Sample is Positive in FAM/Green channel with Control Mix with a CT of 22-30,

If all above are met, then the Test is valid and results can be analyzed further for EGFR mutations.

Above steps for Test validation are summarized in Figure 2.


Fig 2. Test validation flowchart

B) Data analysis for EGFR mutations

- 1) Select the patient’s samples, which are positive in FAM/Green channel with EGFR specific mixes with a CT of 20-40, and document the CTs in the Table 4.
- 2) Calculate ΔCT for the above selected samples through the following equation.

$$\Delta CT = \text{Mutation Mix CT} - \text{Control Mix CT}$$

- 3) Document the ΔCT in the Table 4 too.
- 4) If the ΔCT is within the valid range as mentioned in the Table 4, the sample has the EGFR mutation and Positive.

Above steps for detection of EGFR mutations are summarized in Figure 3.

EGFR Mix	Valid CT	Sample CT	Valid ΔCT	Sample ΔCT	Result
EGFR Ctrl Mix	22-30		-		-
G719X Mix	20-40		≤ 7.4		
19Del Mix	20-40		≤ 7.5		
20Ins Mix	20-40		≤ 6.9		
S768I Mix	20-40		≤ 10.3		
T790M Mix	20-40		≤ 8.2		
L858R Mix	20-40		≤ 10.2		
L861Q Mix	20-40		≤ 10		

Table 4. Valid CT and ΔCT ranges for EGFR Mixes.

As an example, If CT of a sample is 25.3 with Control Mix, 38.4 with G719X Mix, 33.2 with 19Del Mix, 36.8 with 20Ins Mix, 40 with S768I

EGFR RQ (V1.2)

Mix, 36.9 for T790M mix, 40 with L858R Mix and 40 with L861Q Mix, then Δ CT is 13.1 for G719X (38.4-25.3), 7.0 for 19Del (32.3-25.3), 11.5 for 20Ins and 14.7 for S768I, 11.6 for T790M, 14.7 for L858R and 14.7 for L861Q. (Table. 5). According to table 5, delta CT for 19Del is in valid range, and patient is Positive for 19Del.

EGFR Mix	Valid CT	Sample CT	Valid Δ CT	Sample Δ CT	Result
EGFR Ctrl Mix	22-30	25.3	-	-	Valid
G719X Mix	20-40	38.4	≤ 7.4	13.1	Neg
19Del Mix	20-40	32.3	≤ 7.5	7.0	Positive
20Ins Mix	20-40	36.8	≤ 6.9	11.5	Neg
S768I Mix	20-40	40	≤ 10.3	14.7	Neg
T790M Mix	20-40	36.9	≤ 8.2	11.6	Neg
L858R Mix	20-40	40	≤ 10.2	14.7	Neg
L861Q Mix	20-40	40	≤ 10	14.7	Neg

Table 5. EGFR mutation analysis for a specific sample.

Note: above Δ CTs should be validated each run by using Positive and Negative Control which provided in the kit.

Note that a normal sample may cross-react with some of the mutation specific mixes. However, the calculated Δ CT will always fall outside the valid range. Therefore, a sample with CT of 20-40 in FAM/Green channel for EGFR mutation mixes, will only considered positive only if Δ CT is within the valid range.


Fig 3. Sample analysis flowchart for EGFR mutation detection

Results can also be interpreted as below

- Sample is **Negative** for EGFR mutations, if it is negative in FAM/Green channel with all mixes or has CT above 40.
- Sample is **Negative** for EGFR mutations if it is positive in FAM/Green channel with CT of 20-40 for one or more EGFR mixes, but ΔCT is higher than cutoff range (Table 4).

- Sample is **Positive** for EGFR mutations, if sample has a valid CT and Δ CT is lower than cutoff range (Table 4) in FAM/Green channel for one of the seven EGFR mixes. Briefly,
 - It is Positive for G719X if with G719X Mix has a $CT \leq 40$ and $\Delta CT \leq 7.4$. Note that, this kit doesn't differentiate between G719A, G719S and G719C mutations.
 - It is Positive for 19Del if with 19Del Mix has a $CT \leq 40$ and $\Delta CT \leq 7.5$. Note that this kit detects 24 deletion mutations (Table1) but not distinguish between them.
 - It is Positive for 20Ins if with 20Ins Mix has a $CT \leq 40$ and $\Delta CT \leq 6.9$. Note that this kit detects 3 insertion mutations (Table1) but not distinguish between them.
 - It is Positive for S768I if with S768I Mix has a $CT \leq 40$ and $\Delta CT \leq 10.3$.
 - It is Positive for T790M if with T790M Mix has a $CT \leq 40$ and $\Delta CT \leq 8.2$.
 - It is Positive for L858R if with L858R Mix has a $CT \leq 40$ and $\Delta CT \leq 10.2$.
 - It is Positive for L861Q if with L861Q Mix has a $CT \leq 40$ and $\Delta CT \leq 10$.

Interpretation of results for EGFR mutation test are summarized in Table 6.

Note that, if a sample is Negative with this kit in terms of EGFR mutations, the following conditions should be considered:

- Sample is Negative only for EGFR mutations which was mentioned in Table1 and it could be Positive for other EGFR mutations.

- Sample is Positive for EGFR mutations in quantities below the sensitivity of the kit.

EGFR Mix	Sample CT	Sample Δ CT	Conclusion
G719X Mix	>40	-	Neg for G719A/ G719S/ G719C mutations
	20-40	>7.4	Neg for G719A/ G719S/ G719C mutations
	20-40	\leq 7.4	Pos for G719A/ G719S/ G719C mutations
19Del Mix	>40	-	Neg for Exon 19 Deletions
	20-40	>7.5	Neg for Exon 19 Deletions
	20-40	\leq 7.5	Pos for Exon 19 Deletions
20Ins Mix	>40	-	Neg for Exon 20 Insertions
	20-40	>6.9	Neg for Exon 20 Insertions
	20-40	\leq 6.9	Pos for Exon 20 Insertions
S768I Mix	>40	-	Neg for S768I mutation
	20-40	>10.3	Neg for S768I mutation
	20-40	\leq 10.3	Pos for S768I mutation
T790M Mix	>40	-	Neg for T790M mutation
	20-40	>8.2	Neg for T790M mutation
	20-40	\leq 8.2	Pos for T790M mutation
L858R Mix	>40	-	Neg for L858R mutation
	20-40	>10.2	Neg for L858R mutation
	20-40	\leq 10.2	Pos for L858R mutation
L861Q Mix	>40	-	Neg for L861Q mutation
	20-40	>10	Neg for L861Q mutation
	20-40	\leq 10	Pos for L861Q mutation

Table 6. Interpretation of results for EGFR mutation.

13-Analytical Sensitivity

Analytical sensitivity of this assay is equivalent to the percentage of mutated EGFR DNA that can be identified in the background of wild-type DNA. Sensitivity of the assay is determined by dilution of plasmids containing target sequence in human genomic DNA and is mentioned in Table 7. The sensitivity depends on the DNA quality and quantity. Maximum sensitivity is achieved when a sample with Control Mix has a CT of 22-27 in FAM/Green channel.

Reaction	Ctrl Mix CT: 22-27
G719X Mix	2%
19Del Mix	1%
20Ins Mix	2%
S768I Mix	1%
T790M Mix	8%
L858R Mix	1%
L861Q Mix	1%

Table 7. Sensitivity of the EGFR mutation detection assay.