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
# **KRAS RQ Kit Manual**

**For Real-Time PCR Detection of KRAS Mutations**

For use with Rotor-Gene  
Research use only

NG-WI-ASL-33-301  
Version 3.1  
Spring 2022

## Table of Contents:

1. Introduction.....	2
2. Kit Contents.....	3
3. Storage and Stability .....	3
4. General Precautions.....	3
5. Additionally Required Items.....	4
6. Specimen .....	4
7. DNA Extraction .....	4
8. Protocol: DNA Sample Assessment.....	5
9. Programming of Rotor-Gene .....	6
10. Analysis: Sample Assessment.....	6
11. Protocol: Detection of KRAS mutation.....	8
12. KRAS Mutation Detection Analysis .....	10
13. Analytical Sensitivity .....	17

**KRAS RQ** kit is intended for the detection of mutation in codons 12 and 13 of KRAS in human genomic DNA. This kit is designed for use with Rotor-Gene. This kit is for research use only!

## 1- Introduction

Colorectal cancer (CRC) is among the most prevalent cancers worldwide. Treatment with monoclonal antibodies against Epidermal Growth Factor Receptor (Anti-EGFR) has shown to be effective for CRC patients. However, Anti-EGFR would be ineffective in presence of some KRAS mutations including mutations in codons 12 and 13. This is the same for treating non small cell lung cancer (NSCLC) with Anti-EGFR. Therefore, determining the KRAS mutation status is essential for these patients. It should be noted that, 95% and 88% of KRAS mutations in CRC and NSCLC patients respectively are reported in codons 12 and 13 including G12A, G12C, G12D, G12R, G12S, G12V and G13D.

KRAS RQ Kit provides ready to use mixes for detection of these 7 mutations based on Real-Time PCR technology. In this method application of fluorescent probes allows detection of amplified product. Analysis of fluorescent kinetics also leads to detection of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product. This kit also takes advantage of Internal control for preventing false negative results due to failure in DNA extraction or PCR reaction. This kit is designed to be used with Rotor-Gene.

## 2- Kit Contents

The kit contains a manual or a quick guide and following reagents:

Label	Content	Quantity	Volume
Ctrl Mix	PCR Mix for quality control	2	480 µl
G12A Mix	PCR Mix to check G12A mutation	1	480 µl
G12C Mix	PCR Mix to check G12C mutation	1	480 µl
G12D Mix	PCR Mix to check G12D mutation	1	480 µl
G12R Mix	PCR Mix to check G12R mutation	1	480 µl
G12S Mix	PCR Mix to check G12S mutation	1	480 µl
G12V Mix	PCR Mix to check G12V mutation	1	480 µl
G13D Mix	PCR Mix to check G13D mutation	1	480 µl
KRAS Pos	Positive control	1	250 µl
KRAS Neg	Negative control	1	250 µl
Water	PCR Grade Water	1	200 µl

## 3- Storage and Stability

The kit components should be shipped and stored at  $-20^{\circ}\text{C}$  and are stable until the expiry date mentioned. Avoid repeated freeze-thaws especially for PCR Mix more than few times to prevent reduced sensitivity.

## 4- General Precaution

In order to prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation

where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cold blocks instead.

## 5- Additionally Required Items

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Vortex Mixer
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit and required equipments/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

## 6- Specimen

Sections of paraffin block are the proper sample for patients with colorectal cancer and for NSCLC patients, core needle biopsy or aspiration fine needle sample may be used.

## 7- DNA Extraction

DNA extraction can be performed with different kits from various manufacturers. We recommend using:

- QIAamp Fast DNA Tissue Kit (Cat. no. 51404, Qiagen GmbH, Hilden, Germany)
- Roche; High Pure FFPE DNA Isolation Kit (Cat. No. 06 650 767 001, Roche Diagnostics GmbH, Mannheim, Germany)

Note that, the minimum volume of DNA required for the test is 50 microliters. We recommend extracting 100 microliters or higher volume of DNA, enough to repeat the test if necessary.

## 8- Protocol: DNA Sample Assessment

Before examining a sample for KRAS mutations, quality of DNA should be assessed. If the results are within the desired range, then the second test for detecting KRAS mutations will be performed. To qualify DNA extraction, follow steps below.

**Note! Check the volume of extracted DNA sample and make sure it is more than 50 microliters and preferably about 100 microliters. Volumes less than 50 microliters, are not enough to proceed.**

First, thaw the **Ctrl Mix** on ice completely, mix by flickering followed by a quick spin and store on crushed ice after. Place the required number of microtubes on cold block including one for each sample, plus two for positive and negative controls.

**Pipette 20µl of Ctrl Mix to each microtube. Continue by adding 5µl of DNA sample, Positive control or Negative control and water to each tube.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

## 9- Programming Rotor-Gene

*Before you start the machine, make sure the locking ring is attached to the rotor!*

Open the CD provided in the kit and double click on “KRAS 0.2” or on “KRAS 0.1” depending to tubes used. Click on Start button (Green button on the top menu). On the pop up window click start again and save the run file.

You can also use the following table:

Step	Temperature and time	Cycles
1	<b>95°C x 10 min</b>	1
2	<b>95°C x 20 sec</b>	45
	<b>60°C x 60 sec</b>	

Fluorescence should be collected at 60°C for FAM and VIC dyes. Mix contain ROX with final concentration of 300nM in reaction.

## 10-Analysis: Sample Assessment

Data Analysis is performed in two steps. First, the test validity is assessed and then DNA sample analysis. In the first step, results of positive and negative controls will be evaluated.

Perform quantitative analysis for both **DNA sample (Green channel)** and **Internal control (Yellow channel)**. Briefly, click on analysis menu and then, under Quantitation tab double click on “Cycling A. Green”. Close the pop up window and manually set threshold at 0.1. Repeat above for Yellow channel and set the threshold at 0.1.

**Note that, a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample can not be considered as positive (even CT present) result is not reliable.**

**A) Test is valid only if:**

- 1) Water sample is **Negative** in Green channel and it is **Positive** in Yellow channel with CT of 28-31.
- 2) Positive control/Negative control is **Positive** in Green channel with CT of 26-29. And it is **Positive** in the Yellow channel with CT of 28-31.

Expected results for a valid test are summarized in Table1.

Sample	Yellow	Green
NTC	Pos CT: 28-31	Neg
Pos/Neg Ctrl	Pos CT: 28-31	Pos CT: 26-29

Table1. Expected results for Sample Assessment Test validity.

**If above are met, then the Test results are valid and may proceed to DNA sample assessment.**

**B) DNA sample Quality Analysis:**

**DNA sample is qualified only if** sample is positive in Green channel with CT of 22-30 and also positive in Yellow channel with CT of 28-31. This sample can be further examined for KRAS mutations.

**Sample should be diluted with water**, if a sample in is positive in Green channel with CT of less than 22 and also positive in Yellow channel with CT of 28-31. Dilute the sample to reach CT of 22-27. By 2X dilution of sample, CT is increased one cycle. After diluting the sample, proceed to the KRAS mutation test.

**Result is invalid**, if a sample in Green channel has CT above 30. In this case, repeating DNA extraction is recommended.

DNA sample Quality Analysis is summarized in Table 2.

Green	Yellow	Result
+ CT: 22-30	+ CT: 28-31	Valid
+ CT<22	+ CT: 28-31	Sample dilution
- CT>30	CT<28 CT>31	Invalid

Table2. DNA Quality Analysis

## 11-Protocol: Detection of KRAS mutations

To detect KRAS mutations, each sample should be examined with eight mixes, one for DNA quantitation and one for each of seven mutations, each in a separate microtube. Therefore, to examine only one sample, 32 microtubes are required; 8 for the sample, 8 for the Positive control, 8 for the Negative control and 8 for water (NTC). Respectively for each extra samples 8 microtubes will be added. Therefore for 2 samples, 40 microtubes and for 3 samples 48 microtubes are required. Figure1 shows tube setup for 3 samples. Place required number of tubes on cold block organized in series of 8 each.

Pipette 20 $\mu$ l of Ctrl Mix to each tube in the first series.

Pipette 20 $\mu$ l of G12A Mix to each tube in the 2nd series.

Pipette 20 $\mu$ l of G12C Mix to each tube in the 3rd series.

Pipette 20 $\mu$ l of G12D Mix to each tube in the 4th series.

Pipette 20 $\mu$ l of G12R Mix to each tube in the 5th series.

Pipette 20 $\mu$ l of G12S Mix to each tube in the 6th series.

Pipette 20 $\mu$ l of G12V Mix to each tube in the 7th series.

Pipette 20 $\mu$ l of G13D Mix to each tube in the 8th series.

Continue by adding 5 $\mu$ l of extracted DNA, Positive control, Negative Control and water to each tube. Consider the first, second and third tube in each series for positive and Negative control and water. Then next would be for samples.

See Fig 1 for details.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine and program it according to the section 9.

*Note: if using Rotor-Gene attach the locking ring too.*


Fig1. Tubes setup for Mixes and samples.

## 12-KRAS Mutation Detection Analysis

Before analyzing results of KRAS mutations, test validity should be verified. To do so, results of positive and negative controls and samples with Control Mix will be analyzed. If the test is valid, then may proceed to KRAS mutation analysis.

Briefly, click on analysis menu and then, under Quantitation tab double click on “Cycling A. Green”. Close the pop up window and manually set threshold at 0.1. Repeat above for Yellow channel and set the threshold at 0.1.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

**A) Test is valid only if:**

- 1) Water sample is Negative in Green channel with all Mixes, and
- 2) Water sample is Positive in Yellow channel with Control Mix with a CT of 28-31, and
- 3) Positive control is Positive in Green channel with Control Mix with a CT of 26-29, and
- 4) Pos control is Positive in Yellow channel with Control Mix and CT of 28-31, and
- 5) Positive control is Positive in Green channel **with all 7 KRAS mutation Mixes** with a CT of 26-40, and
- 6) Sample is Positive in Green channel with Control Mix with a CT of 22-30 and also Positive in Yellow channel with CT of 28-31.

Above steps for Test validation are summarized in Figure 1.

**If all above are met, then the Test is valid and results can be analyzed further for KRAS mutations.**


Fig 1. Test validation flowchart

## B) Data analysis for KRAS mutations

- 1) Select the patient's samples which are positive in Green channel with KRAS specific mixes with a CT of 20-40, and document the CTs in the Table 3.
- 2) Calculate  $\Delta CT$  for the above selected samples through the following equation.

$$\Delta CT = \text{Mutation Mix CT} - \text{Control Mix CT}$$

- 3) Document the  $\Delta CT$  in the Table 3 too.
- 4) If the  $\Delta CT$  is within the valid range as mentioned in the Table 3, the sample has the KRAS mutation.
- 5) If a sample is positive for two or more KRAS mutations, it is considered positive for the mutation with the lowest  $\Delta CT$  and negative for the other mutations.

Above steps for detection of KRAS mutations are summarized in Figure 2.

KRAS Mix	Valid CT	Sample CT	Valid $\Delta CT$	Sample $\Delta CT$	Result
Ctrl Mix	22-30		NA		
G12A Mix	20-40		$\leq 12$		
G12C Mix	20-40		$\leq 13$		
G12D Mix	20-40		$\leq 8$		
G12R Mix	20-40		$\leq 7$		
G12S Mix	20-40		$\leq 8$		
G12V Mix	20-40		$\leq 3$		
G13D Mix	20-40		$\leq 6$		

Table 3. Valid CT and  $\Delta CT$  ranges for KRAS Mixes.

As an example, If CT of a sample is 26.7 with Control Mix, and 32.4 with G12C Mix, 30.2 with G12R Mix, 34.7 with G12V Mix, 40 with G13D Mix and negative for other mixes, then  $\Delta CT$  is 5.7 for G12C (32.4 - 26.7), 3.5 for G12R (30.2-26.7), 8 for G12V and 13.3 for G13D (Table. 4). So, results of G12C and G12R are in valid range. In this case, patient is Positive for G12R since  $\Delta CT$  is lower for this mutation.

KRAS Mix	Valid CT	Sample CT	Valid $\Delta CT$	Sample $\Delta CT$	Result
Ctrl Mix	22-30	26.7	-	-	Neg
G12A Mix	20-40	Neg	$\leq 12$	-	Neg
G12C Mix	20-40	32.4	$\leq 13$	<b>5.7</b>	<b>Positive</b>
G12D Mix	20-40	Neg	$\leq 8$	-	Neg
G12R Mix	20-40	30.2	$\leq 7$	<b>3.5</b>	<b>Positive</b>
G12S Mix	20-40	Neg	$\leq 8$	-	Neg
G12V Mix	20-40	34.7	$\leq 3$	8	Neg
G13D Mix	20-40	40	$\leq 6$	13.3	Neg

Table 4. KRAS mutation analysis for a specific sample.

Note that a normal sample may cross react with some of the mutation mixes. However, the calculated  $\Delta CT$  will always fall outside the valid range. Therefore, a sample with CT of 20-40 in green channel for KRAS mutation mixes, will only considered positive only if  $\Delta CT$  is within the valid range.


Fig 2. Sample analysis flowchart for KRAS mutation detection

Results can also be interpreted as below;

- Sample is **Negative** for KRAS mutations, if it is negative in Green channel with all mixes or has CT above 40.
- Sample is **Negative** for KRAS mutations if it is positive in Green channel with CT of 20-40 for one or more KRAS mixes, and a  $\Delta CT$  above valid range (which mentioned in Table 3).

- Sample is **Positive** for KRAS mutations, if sample has a valid CT and  $\Delta$ CT (Table 3) in Green channel for one of the seven KRAS mixes. Briefly,
  - It is Positive for G12A if G12A Mix  $CT \leq 40$  and  $\Delta CT \leq 12$ .
  - It is Positive for G12C if G12C Mix  $CT \leq 40$  and  $\Delta CT \leq 13$ .
  - It is Positive for G12D if G12D Mix  $CT \leq 40$  and  $\Delta CT \leq 8$ .
  - It is Positive for G12R if G12R Mix  $CT \leq 40$  and  $\Delta CT \leq 7$ .
  - It is Positive for G12S if G12S Mix  $CT \leq 40$  and  $\Delta CT \leq 8$ .
  - It is Positive for G12V if G12V Mix  $CT \leq 40$  and  $\Delta CT \leq 3$ .
  - It is Positive for G13D if G13D Mix  $CT \leq 40$  and  $\Delta CT \leq 6$ .
- If a sample has a valid CT and  $\Delta$ CT (Table 3) for more than one of the KRAS Mix, sample is Positive for a KRAS mutation which has a lowest  $\Delta$ CT.

Interpretation of results for KRAS mutation test are summarized in Table 5.

Note that, if a sample is Negative with this kit in terms of KRAS mutations, the following conditions should be considered:

- Sample is Negative only for KRAS mutations which was mentioned.
- Sample is Positive for KRAS mutations in quantities below the sensitivity of the kit.
- Sample could be Positive for other KRAS mutations.

KRAS Mix	Sample CT	Sample $\Delta$ CT	Conclusion
G12A Mix	>40	-	Neg for G12A mutation
	20-40	>12	Neg for G12A mutation
	20-40	$\leq$ 12	Pos for G12A mutation
G12C Mix	>40	-	Neg for G12C mutation
	20-40	>13	Neg for G12C mutation
	20-40	$\leq$ 13	Pos for G12C mutation
G12D Mix	>40	-	Neg for G12D mutation
	20-40	>8	Neg for G12D mutation
	20-40	$\leq$ 8	Pos for G12D mutation
G12R Mix	>40	-	Neg for G12R mutation
	20-40	>7	Neg for G12R mutation
	20-40	$\leq$ 7	Pos for G12R mutation
G12S Mix	>40	-	Neg for G12S mutation
	20-40	>8	Neg for G12S mutation
	20-40	$\leq$ 8	Pos for G12S mutation
G12V Mix	>40	-	Neg for G12V mutation
	20-40	>3	Neg for G12V mutation
	20-40	$\leq$ 3	Pos for G12V mutation
G13D Mix	>40	-	Neg for G13D mutation
	20-40	>6	Neg for G13D mutation
	20-40	$\leq$ 6	Pos for G13D mutation

Table 5. Interpretation of results for KRAS mutation.

### 13-Analytical Sensitivity

Analytical sensitivity of this assay is equivalent to the percentage of mutated KRAS DNA that can be identified in the background of wild-type DNA and is mentioned in Table 6. The sensitivity depends on the DNA quantity in a sample. Maximum sensitivity of is obtained when a sample with Control Mix has a CT of 22-27 in

Green channel. If CT range is 28-30, the sensitivity will decrease for some mutations.

<b>Reaction</b>	<b>Ctrl Mix CT: 22-27</b>
G12A	1%
G12C	4%
G12D	1%
G12R	2%
G12S	1%
G12V	4%
G13D	4%

Table 6. Sensitivity of the KRAS mutation detection assay.