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Factor XIII

RQ Kit Manual

For Real-Time PCR Detection of
Factor XIII (Val34Leu)

For use with Rotor-Gene 6000/Q
Or StepOne/StepOnePlus

Research use only

NG-WI-ASL-19-100
Version 1.0
March 2017


Table of Contents:

1. Introduction	2
2. Kit Contents.....	3
3. Storage and Stability	3
4. General Precautions	3
5. Additionally Required Materials	4
6. Specimen, Storage and Transport.....	4
7. Interfering Substances	5
8. DNA Isolation	5
9. Internal Control.....	5
10. PCR Protocol	5
11. Programming RotorGene.....	6
12. Programming StepOne.....	6
13. Programming Other Machines	7
14. Data Analysis: RotorGene	7
15. Data Analysis: StepOne	9

FXIII RQ

Factor XIII RQ kit is intended for the detection of Val34Leu polymorphism of Factor XIII in genomic DNA extracted from human cells. This kit is designed for use with RotorGene 6000/Q or StepOne/StepOnePlus instrument. This kit is for research use only!

1. Introduction

Factor XIII, first called fibrin stabilizing factor, was discovered in 1948. It is involved in the last stage of coagulation and participates in clot preservation. It also participates in wound repair and healing, maintaining pregnancy and acts as proangiogenesis. Among the few polymorphisms known for this factor, the Val34Leu has been studied more extensively as it increases the activity of this factor and leads to change in fibrin clots. Studies have shown that this polymorphism may reduce the risk of myocardial infarction, brain infarction and venous thromboembolism (VTE). However, it also increases the risk of intracerebral hemorrhage and pulmonary embolism.

Factor XIII RQ kit provides a ready-to-use Real-Time PCR system for detection of Val34Leu polymorphism of Factor XIII with RotorGene 6000/Q or StepOne/StepOnePlus instrument. This method applies fluorescent dye labeled probes allowing detection of amplified product. Analysis of fluorescent kinetics also leads to assessment of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

FXIII RQ

2. Kit Contents

The kit contains a manual, a CD with RotorGene and StepOne templates and following reagents:

Label	Content	Quantity
F13Q Mix	PCR Mix*	500 μ l
F13 MM Ctrl	Homozygous Positive Control	100 μ l
F13 WM Ctrl	Heterozygous Positive Control	100 μ l
F13 WW Ctrl	Negative Control	100 μ l
Water	PCR Grade Water	200 μ l

*1, 2 or 4 tubes for 24, 48 or 96 reaction kits, respectively.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws especially for PCR Mix more than few times to prevent reduced sensitivity.

4. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.

FXIII RQ

- Thaw kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Keep PCRMix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cooling blocks instead.

5. Additionally Required Materials

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex mixer
- Adjustable pipetters
- DNA extraction kit
- Nuclease free filtered tips
- Tubes or strips and strip caps
- Disposable powder-free gloves
- Cold block

6. Specimen, Storage and Transport

Whole blood (0.5ml) is the preferred sample. Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at $+4^{\circ}\text{C}$ (stable for few days). For longer terms, sample should be aliquoted and stored at -20°C which is stable for few months.

FXIII RQ

7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinised patients must not be used as well.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

8. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using following:

- High Pure PCR Template Preparation Kit (Cat. no. 11796828001, Roche Applied Science, Mannheim, Germany).
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany).

9. Internal Control

This kit requires no internal control since each person carries either wild type or mutant or both alleles. Therefore, any sample is always positive for at least one of them. If a sample is negative for both alleles, then the test should be repeated.

10. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cooled PCR block. Consider one tube for each sample plus one for each control and one for the NTC.

FXIII RQ

Pipette 20ul of F13Q Mix directly to each tube followed by adding 5ul of controls or isolated DNA.

Cap the tubes/Strips and visually inspect to make sure all are capped securely. If using StepOne/StepOnePlus instrument, spin tubes/strips briefly. Place tubes in the machine. When using RotorGene attach the locking ring.

11. Programming RotorGene

- Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the CD provided in the kit and double click on F13Q0.2 (if 0.2ml tubes are used) or on F13Q Strip (if 0.1ml tubes are used) to open the program. Click on Start button (Green button on the top menu). On the pop up window click start again and save the run file.

Edit sample names.

12. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on "Plate Setup". Controls and few samples are defined. You may change plate set up using right click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished, click on Start Run and save the experiment on desired location. Instrument will start shortly.

FXIII RQ

13. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 15 sec	40
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. PCR Mix contains ROX with the final concentration of 300nM in the reaction.

14. Data Analysis: RotorGene

Analyze data according to manufacturer recommendations. Briefly, click on analysis menu select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both Green and Yellow channels and then click on “Show”.

The “Scatter Graph Window” is now shown in the lower left window. Each dot represents one sample. The vertical axis shows the Green fluorescence and horizontal axis is for the Yellow fluorescence. Note that M allele (Leu 34) is detected in Green channel and W (Val34) allele in Yellow channel.

Therefore, samples are located in four regions of the graph. MM samples with high Green and low Yellow fluorescence gather in upper left; WW samples with low Green fluorescence and high Yellow gather in lower right; WM samples with high Green and high Yellow fluorescence gather in upper right; and finally NTC or

FXIII RQ

control with no DNA gather in lower left with low Green and low Yellow fluorescence.

In order to determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then drag a rectangle around each control. Each rectangle represents one of the above mentioned regions. Label MM genotype region as “Mutant”, WM as “Heterozygous”, WW as “Wild Type” and Negative Control as “None”(Figure 1).


Figure 1: Defining genotypes of different controls of the test on Rotor-Gene instrument.

Now turn on all the samples and the results are displayed in the tab “Scatter Analysis Results” (Figure 2).


Figure 2: Typical scatter graph for controls and samples on Rotor-Gene instrument.

Note! If location of controls are not similar to the typical graph represented here; or if any of above regions overlap in a graph, results are not valid and test should be repeated. Also if a sample is located in between regions, it should be re-examined.

15. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and select “Allelic Discrimination” tab. Genotype of the samples are automatically determined by the software based on comparison of their fluorescence with the controls.

Figure 3 represents the typical graph. Each dot represents one sample. All samples are color coded based on their genotype as blue for MM, Red for WW, Green for WM, and Black box for control with no DNA or NTC. Unknown samples will be marked with black "×".

Note that the vertical axis shows the FAM fluorescence and horizontal axis is for the VIC fluorescence. Accordingly M allele is detected in FAM channel and W allele in VIC channel. Therefore, samples are located in four regions of the graph. MM samples with high FAM and low VIC fluorescence gather in upper left; WW samples with low FAM fluorescence and high VIC gather in lower right; WM samples with high FAM and high VIC fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low FAM and low VIC fluorescence. Genotypes of samples are automatically determined by the software based on comparison with the controls (Figure 3).


Figure 3: Typical graph for controls and samples on StepOne instrument.

Note! If location of controls are not similar to the typical graph represented here; or if controls are located too close to each other in a graph; and if samples lie in between controls and genotypes cannot be called, results are invalid and test should be repeated.